


# ProTime<sup>®</sup>

Microcoagulation System

**NEW**

**GET  
CONNECTED!  
EMR READY**



**CLIA WAIVED**

## Advanced Technology for Improved Patient Care

- Accepts and stores patient and/or operator ID for reduced transcription errors
- Auto-send results quickly and easily to printer, computer or LIS
- Only system offering the flexibility of both on-board and optional external quality controls
- Thromboplastin reagent of ISI = 1.0 (AHA, CAP, and WHO recommended)
- Off-meter sampling provides safety against disease transmission
- Barcoded cuvette - no coding necessary
- Longer open-pouch stability, up to 16 hours
- Available for patient self-testing




THE POINT OF CARE

## ORDERING INFORMATION

Catalog	Description	Packaging
PROTIMEPRO	ProTime professional instrument, AC power module, AC power cord, operator's manual, training CD, 3-year warranty	1 each
PRO3-25	ProTime disposable 3-channel cuvettes (low blood volume) consists of: 25 cuvettes 25 Tenderlett Plus LV Package insert	1 box; 10 boxes/case
PRO3-50	ProTime disposable 3-channel cuvettes (low blood volume) consists of: 50 cuvettes 50 Tenderlett Plus LV Package insert	1 box each: cuvettes Tenderletts
PRO5-25	ProTime disposable 5-channel cuvettes consists of: 25 cuvettes 25 Tenderlett Plus Package insert	1 box; 10 boxes/case
DCPRO-N	Liquid controls for ProTime - Normal (optional)	1 box of 15 ampules
DCPRO-A	Liquid controls for ProTime - Abnormal (optional)	1 box of 15 ampules
LBLKIT	Printer, label maker	1 each
LBLPAPER	Paper for LBLKIT	1 each
TL25PRO3	Tenderlett Plus LV finger incision device	1 case; 18 bags/25 TL+ LV
TL25PRO5	Tenderlett Plus finger incision device	1 case; 18 bags/25 TL+
PROBAG	ProTime instrument carrying bag	1 each
PROCABLE	ProTime computer cable	1 each
JavaLin 300*	Clinical interface	

\* Contact Dawning Technologies, Inc. at 1-800-332-0499 or [www.dawning.com](http://www.dawning.com)


### THREE EASY STEPS



**1. Collect Sample**



**2. Attach Tenderlett® Plus**



**3. Read Results**

**1-800-631-5945**

**[www.protimesystem.com](http://www.protimesystem.com)**

ProTime and Tenderlett are registered trademarks of ITC.



8 Olsen Avenue • Edison, NJ 08820 USA  
tel: 732.548.5700 • fax: 732.248.1928  
[www.itcmed.com](http://www.itcmed.com)

a subsidiary of Thoratec Corporation



ISO 9001 Certified

CPT Code: 85610QW  
E&M Code: 99211