






DEALER QUICK REFERENCE GUIDE

PRODUCT	DESCRIPTION	ORDER NUMBER	CONTENTS
 <p>CLIA WAIVED</p>	<p>ProTime professional instrument</p> <p>Features:</p> <ul style="list-style-type: none"> • Accepts and stores patient and/or operator ID for reduced transcription errors • Auto-send results quickly and easily to printer, computer or LIS • Only system offering the flexibility of both on-board and optional external quality controls <p>CPT Code 85610QW, average Medicare reimbursement rate \$5.49 CPT Code (E&M) 99211, average Medicare reimbursement rate \$21.60</p>	PROTIMEPRO	<p>Instrument</p> <p>Training CD</p> <p>Charger</p> <p>Operator's Manual</p> <p>3 year warranty</p> <p>Carrying bag</p>
	<p>ProTime3 and 5-channel cuvettes</p> <p>Both perform on-board Level 1 and Level 2 controls with each patient test. ProTime 5-channel offers triplicate testing. ProTime 3-channel offers the same accuracy and reliability of the 5-channel using one, large drop of blood.</p> <ul style="list-style-type: none"> • Bar-coded cuvette - no coding necessary • Thromboplastin reagent of ISI = 1.0 (AHA, CAP, and WHO recommended) <p>STABILITY: Cuvettes are ready to use. Unopened cuvettes stored at (2-8°C, 35-46°F) are stable until date printed on pouch, or for 60 days at room temperature. Open pouch stability is 16 hours.</p>	<p>PRO3-25</p> <p>PRO5-25</p> <p>PRO3-50</p>	<p>ProTime3 (25 cuvettes, 25 Tenderlett Plus LV)</p> <p>ProTime (25 cuvettes, 25 Tenderlett Plus)</p> <p>ProTime3 (50 cuvettes, 50 Tenderlett Plus LV)</p>
	<p>Tenderlett® Plus LV incision device - Included with cuvettes</p> <p>Two-in-one incision and sample collection device included with the ProTime 3-channel cuvettes. Tenderlett Plus is supplied with the ProTime 5-channel cuvettes. Both sizes are also sold separately.</p>	<p>TL25PRO3</p> <p>TL25PRO5</p>	<p>Case = 18 bags (25 devices per bag)</p>
	<p>ProTime printer, label maker</p> <p>Connects directly to the ProTime instrument. Results are printed on a label that can be easily applied to the patient's chart. Additional label paper available separately.</p>	<p>LBLKIT</p> <p>LBLPAPER</p>	<p>1 unit</p> <p>2 per package</p>
	<p>ProTime liquid controls</p> <p>Optional. Generally not required with on-board controls for CLIA-waived testing.</p>	<p>DCPRO-N</p> <p>DCPRO-A</p>	<p>Normal (1 box of 15 ampules)</p> <p>Abnormal (1 box of 15 ampules)</p>

ProTime® and Tenderlett® are registered trademarks of ITC.